



## **Titration of HCL with NaOH**

Kimmy Pathanasap (5861052)

Supichaya Nobnom (5861155)

Pichayathida Siriwechdaruk (5861105)

Pnotporn Jantarakolica (5861115)

Passawas Satitjindawong (5861095)

Advanced Chemistry 1102

Ms.Patrachorn Sanguansat


Mahidol University International Demonstration School

Semester 2 Academic year 2016-17

## Discussion

### Questions:

1. How would it affect your results if you used a beaker with residual water in it to measure out your standardized sodium hydroxide solution?
  - Our experiment would not be accurate because there is initially some leftover water inside, which would affect the pH of the acid inside the erlenmeyer flask.
2. How would it affect your results if you used a wet Erlenmeyer flask instead of a dry one when transferring your acid solution from the volumetric pipette?
  - With a wet Erlenmeyer flask, the result will not be accurate. It might be contaminated by any type of liquid such as basic or acidic. Therefore, an error will occur since it is not 100% acid.
3. How do you tell if you have exceeded the equivalence point in your titration?
  - We have to find an indicator which will be able to tell us whether the solution is neutralized into pH 7. The indicator might have different color for each pH. So, it will be able to tell us the pH of that particular solution. For example, Bromthymol Blue which its color will be turned to blue when the solution is acidic, green if it solution is neutral, and blue if the solution is basic.
4. Vinegar is a solution of acetic acid ( $\text{CH}_3\text{COOH}$ ) in water. For quality control purposes, It can be titrated using sodium hydroxide to assure a specific % composition. If 25.00 mL of acetic acid is titrated with 9.08 mL of a standardized 2.293 M sodium hydroxide solution, what is the molarity of the vinegar?



$$\text{mol NaOH}(aq) = \text{concentration of NaOH} \times \text{volume of NaOH} = 0.00908 \times 2.293 = 0.02082044 \text{ mol}$$


$$0.02082044 \text{ mol NaOH} = 0.02082044 \text{ mol CH}_3\text{COOH}$$


$$\text{CH}_3\text{COOH concentration} = \frac{0.02082044}{0.025} = 0.8328176 \text{ mol/L}$$

- The molarity of the vinegar is 0.83.

### Why result turn out in different shades of pink?

- The different shades of pink color is based on the amount of NaOH released into the beaker. If there is an increase in NaOH, then the color of pink will become darker. However, if the amount of NaOH is decreased, then the color of pink will be lighter.



**Calculation:**

$$\text{Molarity} = \frac{\text{mol}}{\text{volume (L)}}$$

$$\text{mol base} = \text{molarity} \times \text{volume} = 0.1 \times 0.009 = 0.0009 \text{ mol acid}$$

From equation: Acid to Base mole ratio = 1:1

Thus, mol acid = mol base = 0.0009 mol

$$\text{Molarity of base} = \frac{\text{mol}}{\text{volume (L)}} = \frac{0.00009}{0.0010} = 0.009 \text{ mol/L}$$

$$\text{Average acid concentration} = \frac{\Sigma \text{concentration of acid}}{\Sigma \text{number of trials}} = \frac{0.0087+0.009+0.0089+0.0091+0.0094}{4} = 0.0074 \text{ mol/L}$$

|                                    | Trial 1 | Trial 2 | Trial 3 | Trial 4 | Trial 5 |
|------------------------------------|---------|---------|---------|---------|---------|
| Volume of base (L)                 | 0.0087  | 0.009   | 0.0089  | 0.0091  | 0.0094  |
| Concentration of base (mol/L)      | 0.1     | 0.1     | 0.1     | 0.1     | 0.1     |
| Mole of base (mol)                 | 0.00087 | 0.0009  | 0.00089 | 0.00091 | 0.00094 |
| Acid to Base (mol)                 | 1:1     | 1:1     | 1:1     | 1:1     | 1:1     |
| Moles of acid (mol)                | 0.00087 | 0.0009  | 0.00089 | 0.00091 | 0.00094 |
| Volume of acid (L)                 | 0.01    | 0.01    | 0.01    | 0.01    | 0.01    |
| Acid concentration (mol/L)         | 0.087   | 0.09    | 0.0089  | 0.0091  | 0.0094  |
| Average acid concentration (mol/L) | 0.0074  |         |         |         |         |

**Conclusion**

According to the experiment, 0.009 litres of 0.1M sodium hydroxide can neutralize 0.01 litres of 0.09M hydrochloric acid effectively.